



Dog CPR Flowchart

**Check for breathing
for 10 seconds**

- Open airway
- Check for normal breathing for no more than 10 seconds

Call for help if available

- Call for help and get someone to call your vet if possible

Deliver 5 initial rescue breaths

- Open airway, pull the tongue to one side and seal the mouth. Deliver 5 rescue breaths until you see the chest rise by blowing in the nose

Deliver 30 chest compressions

- Deliver 30 compressions in centre of chest at a rate of 100-120 compressions per minute to a depth of at least one third the depth of the chest with one hand

Deliver 2 rescue breaths

- Then open the airway and deliver 2 breaths until you see the chest rise

Continue CPR 30:2

- Repeat at 30 compressions and 2 breaths

**Continue until they start to
breathe**

- Continue until they start to breathe, you are too tired to continue, someone else takes over or the vet arrives.

ProTrainings Europe Ltd

01206 805359

www.protrainings.eu

www.propetfirstaid.co.uk